

# IDI MDS™

## Mobile Display & Support System

### Ergonomic design

Arm height adjustment of 36" means comfortable viewing in both seated and standing positions

#### Easily position:

- Monitors
- Lights
- Radiation Shields



### Mobile imaging without compromise

- Single and dual-arm configurations available
- Each arm provides 36" of fully counter-balanced vertical travel\*
- Low-force manual motion translates into easy positioning and patient safety
- Monitors & support equipment where you need them and up and out-of-the-way when you're done
- Built-in adjustable equipment racks; image sources and monitors can stay connected
- Easy to move with large lockable double-wheeled casters

## Your imaging system & procedure table is mobile... Shouldn't your support equipment be as well?

The IDI MDS mobile display and support system provides the benefits you're used to in fixed-room support equipment; now available in a mobile platform.

*Make the switch to fully mobile imaging without compromise.*

idi

Image Diagnostics, Inc.

Safety, Simplicity, Quality

\*Patent pending

# IDI MDS

## Specifications:

- Retracted height (chassis) : 64.5" (163.8 cm)
- Extended arm vertical height: 96.0" (243.8 cm)
- Arm(s) radius extended horizontally: ~ 76" (193 cm)
- Arm capacity: 48 lbs. (21.8 kg) each
- Monitor mounts: 100 mm<sup>2</sup> and 100 mm x 200 mm VESA mounts available

## Image Diagnostics

Since 1986, IDI has been a manufacturer of innovative products for surgical, interventional and diagnostic imaging.

Our product focus is in electro-mechanical solutions for fluoroscopic and digital radiographic applications.

IDI is committed to superior product design, value and customer support in everything we do.

### Image Diagnostics, Inc.

310 Authority Drive  
Fitchburg, MA, 01420 USA  
Toll free: 877-304-5434  
Tel: 978-829-0009  
Fax: 978-829-0027  
www.imagediagnostics.com



Monitors shown in vertically extended and compact positions

MDS with IDI ISR table




Image Diagnostics, Inc. reserves the right to update product specifications, without notification, as required.